This 24 hour collection kit MUST be mailed to the lab as soon as possible.





Please check all before shipping:

- \square Placed tubes in kit box
- \square Placed frozen ice pack on tubes in kit box
- Completed questionnaire (online or paper form)
- □ Indicated collection times (online or paper form)
- \square Placed completed questionnaire in kit box
- □ Placed completed billing/PHI form in kit box
- \square Submitted a copy of front & back of insurance card (online or in kit)

Thank you for choosing Physicians Lab! 877-316-8686

© Copyright 2018 Physicians Lab. All rights reserved. PHYL-KIT-FEMALE-24HR-2018-07-17

 \downarrow Tear here and return form below with kit \downarrow

Here's she we here's the sheet

Symptoms: Female Patients

Abdominal pain					Limited flexibility				
Abundant light-colored urine					Longer to recover after exercise				
Aches and pains					Loss of body and pubic hair				
Acne					Loss of scalp hair				
Allergies					Low blood pressure				
Anxious					Low blood sugar				
Bleeding changes					Memory lapse				
Bone Loss					Mood swings				
Bruises Easily					Morning fatigue				
Cold body temperature					Nails breaking or brittle				
Constipation					Nervous				
Coughing up blood					Night sweats				
Coughing up tinged phlegm					Numbness–feet or hands				
Craves salt and salty foods					Oily skin or hair				
Decreased libido	Π		Π		Overheating easily	Π	Π	Π	Π
Decreased muscle size	Ē	Π	Π	Ē	Painful intercourse	Ē	Π	Π	Ē
Decreased stamina	П	Π	П		Pale face	П	П	П	Π
Decreased sweating	п	П	П	F	Poor appetite	Ħ.	П	T	T I
Deep forehead wrinkles	П		П		Poor memory	П	П	П	П
Depressed	Ы	П	Ы	F	Prolonged periods	H.	П	П	П
Diarrhea	Ы	П	Н		Prop yourself up to sleep better	Н	Н	Н	П
Does not feel rested in morning	Ы	H	Ы	F	Rapid aging	H.	Ы	Н	П
Easily fatigued when exercising	Ы	П	Ы		Rapid heartbeat	Н	П	П	П
Elevated triglycerides	н	H	н	F	Sagging cheeks	H	H	H	H
Evening fatigue	Н	Н	Н	F	Sees color less vividly	H	Н	Н	Н
Excessive coughing	H	H	н	F	Sensitivity to chemicals	H	H	H	H
Excessive sweating	Н	Н	Н		Sensitivity to heat	H	Н	H	Н
Eyebrows thinner on the sides	Н	Н	Н	F	Shortness of breath	H	H	H	H
Fibrocystic breasts	H	Η	H	F	Sleep disturbed	H	H	H	H
Fibromyalgia	H	H	н		Slow pulse rate	H	H	H	H
Foggy thinking	H	H	H		Snoring	H	H	H	H
Frequent respiratory infections	H	H	Η	F	Strange dreams	H	H	H	H
Gains weight easily	H	H	H	F	Stress	H	H	H	H
Goiter	H	Н	Н	F	Sugar craving	H	Η	H	H
Hair dry or brittle	H	H	H	F	Swelling of ankles	H	Η	H	H
Headaches	H	H	Η	1	Swelling or puffy eyes, face	H	Η	H	H
Hearing loss	H	H	H		Tearful	H	H	H	H
Heart palpitations	H	H	H	-	Tender breasts	H	H	H	H
High blood pressure	H	H	H		Thinner lips	H	H	H	H
	H	Η	Н	F	Thinning skin	H	H	H	H
High cholesterol Hoarseness	H	H	H	-	Thirsty	H	H	H	H
Hot flashes	H	H	Η	-	- ,	H	H	H	H
Incontinence	H	H	H		Trouble falling asleep if woken Unable to lose fat after diet	H	H	H	H
	H	H	H	-	Unable to lose fat after exercise	H	H	H	H
Increase facial or body hair	H	님	H	-		H	H	H	H
Increased appetite	H	H	H	는	Uses alcohol to help fall asleep	H	H	H	H
	H	님	H	-	Uses meds to help fall asleep	H	H	H	H
Increased number of age spots	H	H	H		Uterine fibroids	H	H	H	H
Increased sex drive	H	H	H		Vaginal dryness	H	H	H	H
Increased urinary urge	H	H	H	F	Warm hands and feet	H	H	H	H
Infertility problems	님		님		Water retention	님	님	님	님
Irritable	H		H		Weight gain—hips	H	H	H	H
Joint pain	님		님		Weight gain–waist	님	님	님	님
Joints swelling	H	H	H	F	Wheezing				
Lightheaded when standing					l				

PLEASE TEAR OFF THIS SHEET AND RETURN IT WITH YOUR TESTING KIT. YOU MUST ALSO RETURN THE BILLING/PHI FORM AND A COPY OF THE FRONT AND BACK OF YOUR INSURANCE CARD TO ENSURE YOUR SAMPLES ARE PROCESSED. Thank you for choosing Physicians Lab as your urinary hormone testing provider.

FEMALE: 24 HOUR COLLECTION INSTRUCTIONS IMPORTANT: Confirm that it is YOUR NAME on all collection tubes. Please call 877-316-8686 if your name is NOT listed on the collection tubes. TAKE ALL MEDICATIONS AS PRESCRIBED! See exceptions* **DO NOT STOP** *If you are taking Oral DHEA, Oral Testosterone, and/or Oral Estrogen **TAKING YOUR** Take 16 hours before collecting and do not take the day of collection. The day after your collections are complete, resume as prescribed. **HORMONES!** TIMING GUIDE - CHOOSING YOUR COLLECTION DATE O If you experience a menstrual cycle: O If you use hormone creams/topical hormones on the genitals: • The optimal time to collect is 7 days before you expect your next period. Apply your hormones as you normally do until the day prior to testing. O If you administer or receive injections/pellets • The day prior to testing - apply your hormone cream to your upper inner Start collecting at midpoint between injection/pellet inserts. arm and rub in well. **O** If you use patches • On the day of testing - wake up and use a "Clean Catch" technique to collect • Start collecting 1-2 days after applying patch. your sample then apply your cream to your upper inner arm. Rub in well and O If you use melatonin, hydrocortisone, and/or pregnenolone: be certain to use gloves or an applicator to apply hormones during testing. • Start collecting at least 16 hours after the last dose. • Continue to collect using a "clean catch" technique for the remainder **O** If none of these apply to you: of the day. You may collect on any day **unless** you are taking hormones. If you are taking Clean Catch Technique - Wipe from front to back with a clean damp cloth. Hold the hormones, collect on days you take hormones. labia apart while collecting urine. This will help prevent from washing the vaginal walls and contaminating the urine with residual cream. **IMPORTANT:** Most rejected samples are a result of exceeding fluid restrictions, consuming caffeine and/or taking diuretic medications. Please drink less than 8 oz. in between each collection and minimize caffeine/diuretics to avoid sample dilution. Dilute samples will require re-collection and may incur an additional shipping fee. All collection tubes must be sent back or samples may be rejected.

If you are unsure about the exceptions and/or timing guide above, please contact our staff at 877-316-8686 for assistance before collecting.

Instructions for Collection

The night before testing:

- Limit your caffeine and diuretic intake for the next 24 hours and stop drinking fluids between 8pm and your first morning collection.
- Completely empty your bladder before bedtime. Do not collect this urine.
- Place the ice pack into the freezer.

On the day of collection:

- Collect urine into a disposable cup (not provided). Transfer urine from cup to color coded tubes as indicated below using plastic pipette (included).
- You will need to limit yourself to 5 urination times that will be collected 4 hours apart. If you are unable to limit the number of urinations to 5, call 877-316-8686 for further instructions. You may also use our text message reminder system to schedule collection times and prompt you when it is time to collect.
- Log your urination times on the page below where indicated. You may also enter your collection times online by logging in to PhysiciansLab.com, or by using our email and text messaging instructions.
- Samples can be kept at room temperature for a maximum of 12 hours, after which they MUST be refrigerated.
- Keep urine samples in your refrigerator until they are ready to be shipped.

Shipping:

• Place the urine samples and frozen ice pack into the box immediately before shipping and return using the pre-paid shipping envelope provided.

Overnight Collection: If you wake in the middle of the night to urinate, collect this sample in the **yellow-top** tube labeled "Overnight Collection." If you do not wake up in the middle of the night, fill both yellow tubes with your waking collection sample.

Collection 1: Collection 2: Collection 3: Collection 4: Collection 5: Collect your first sample 4 hours later, collect your in the **yellow-top** second sample in the third sample in the fourth sample in the **pink/** fifth and final sample in the tube labeled "waking orange-top tube labeled red-top tube labeled lavender-top tube labeled **blue-top** tube labeled collection" upon waking. "4 hours after waking." "8 hours after waking." "12 hours after waking." "16 hours after waking." You may collect on any day of the week, but **Returning Your Kit:** must ship the samples between Mon-Sat ONLY. Place all these items in box: In the event you are not able to ship them next day, Tubes with frozen simply store the samples in the **refrigerator** until you can ship them. ice pack on top Place your kit in the FedEx envelope and bring it to any certified • Billing/PHI and **FedEx** FedEx drop off location, drop box, or local FedEx Office. questionnaire forms Place your kit in the USPS envelope and place it back in your • Copy of your **UNITED STATES** POSTAL SERVICE • mailbox for pick up, or bring it to any USPS drop box or local insurance card I AB (front & back) USPS office. Need more help? Call 877-316-8686 Monday–Friday, 8am–5pm EST Tear here and return form below with kit

EMALE OUESTIONNAIRE: GENERAL INFORMATION

TEMALE QUESTIONNAILE. GENERAL INFORMATION

You may skip this page if: You complete this questionnaire online at PhysiciansLab.com or through our text messaging and email instructions.


Check here if completed online: 🗌

Name:
Height:/ feet/inches Weight: pounds Currently Pregnant? 🗌 No 🗌 Yes Due Date:// mm/dd/yy
Hysterectomy? No Yes Hysterectomy Year: 4 digit year Ovaries Removed? N/A One Both Ovaries Removed Year: 4 digit yr
Menstrual Cycles: 🗌 Regular 🗌 Irregular 🗌 None 🛛 First Day of Last Menses:/

Collection Day: ____/___/ ___ mm/dd/yy

Please log your collection time below.

See instruction page for details on when and how to collect, as well as contact information for our customer support team if you have any questions.



Medication Questions

Currently taking inflammatory medication containing cortisone or other anti-inflammatory steroids? 🗌 No 🗌 Yes 🛛 I am not currently taking any hormones 🗌

Hormone Medications

The medication information collected in the section below is essential for the accuracy of your laboratory results and will be displayed therein. Please check off all applicable information. If you select one or more Delivery Method(s) please indicate how long you have been taking each hormone.

Delivery Method	Testosterone*	Progesterone*	Estradiol*	DHEA*	Pregnenolone*	Estriol*	Adrenal Supplement*
Oral							
Injection							
Pellet							
Cream/Gel							
Troche							
Sublingual							
Patch							
How Long?							
0-3 Months							
3-6 Months							
6-12 Months							
12 +							

*If you are taking a medication that combines any of the hormones listed above, please check each of those individual hormones.





© Copyright 2018 Physicians Lab. All rights reserved. PHYL-KIT-FEMALE-24HR-2018-07-17