to the lab as soon as possible. This collection kit MUST be mailed

State-of-the-art Science. Superior Solutions.



9898-912-778 Thank you for choosing Physicians Lab!

□ Submitted a copy of front & back of insurance card (online or in kit ☐ Placed completed billing/PHI form in kit box □ Placed completed questionnaire in kit box □ Indicated collection times (online or paper form) □ Completed questionnaire (online or paper form) ☐ Placed frozen ice pack on tubes in kit box □ Placed tubes in kit box

Please check all before shipping:



Returning Your Kit:

Place all these items in box:

- Tube with frozen ice pack on top
- Billing/PHI and questionnaire forms
- Copy of your insurance card (front & back)



You may collect on any day of the week, but must ship the samples between Mon-Sat ONLY.

In the event you are not able to ship them next day, simply store the samples in the **refrigerator** until you can ship them.



Place your kit in the FedEx envelope and bring it to any certified FedEx drop off location, drop box, or local FedEx Office.



Place your kit in the USPS envelope and place it back in your mailbox for pick up, or bring it to any USPS drop box or local

Need more help? Call 877-316-8686 Monday-Friday, 8am-5pm EST

	↓ Tear here and retur	n form below with kit 🗼	
Symptoms:	d silve d se del		Herel chillips of the higher the comment
Female Patients	Here's child theyere where		Here they office they are
Abdominal pain		Limited flexibility	
Abundant light-colored urine		Longer to recover after exercise	HHHH
Aches and pains		Loss of body and pubic hair	
Acne		Loss of scalp hair	
Allergies		Low blood pressure	
Anxious		Low blood sugar	HHHH
Bleeding changes		Memory lapse	
Bone Loss		Mood swings	
Bruises Easily		Morning fatigue	
Cold body temperature		Nails breaking or brittle	
Constipation		Nervous	
Coughing up blood	百百百百	Night sweats	百百百百
Coughing up tinged phlegm	пппп	Numbness—feet or hands	
Craves salt and salty foods	T T T T	Oily skin or hair	7 7 7 7 7 T
Decreased libido	пппп	Overheating easily	
Decreased muscle size		Painful intercourse	
Decreased stamina		Pale face	
Decreased sweating		Poor appetite	
Deep forehead wrinkles		Poor memory	
Depressed		Prolonged periods	
Diarrhea		Prop yourself up to sleep better	
Does not feel rested in morning		Rapid aging	
Easily fatigued when exercising		Rapid heartbeat	
Elevated triglycerides		Sagging cheeks	
Evening fatigue		Sees color less vividly	
Excessive coughing		Sensitivity to chemicals	
Excessive sweating		Sensitivity to heat	
Eyebrows thinner on the sides		Shortness of breath	
Fibrocystic breasts		Sleep disturbed	
Fibromyalgia		Slow pulse rate	
Foggy thinking		Snoring	
Frequent respiratory infections		Strange dreams	
Gains weight easily		Stress	
Goiter		Sugar craving	
Hair dry or brittle		Swelling of ankles	
Headaches		Swelling or puffy eyes, face	
Hearing loss		Tearful	
Heart palpitations		Tender breasts	
High blood pressure		Thinner lips	
High cholesterol		Thinning skin	
Hoarseness		Thirsty	
Hot flashes		Trouble falling asleep if woken	
Incontinence		Unable to lose fat after diet	
Increase facial or body hair		Unable to lose fat after exercise	
Increased appetite		Uses alcohol to help fall asleep	
Increased fat		Uses meds to help fall asleep	
Increased number of age spots		Uterine fibroids	
Increased sex drive		Vaginal dryness	
Increased urinary urge		Warm hands and feet	
Infertility problems		Water retention	
Irritable		Weight gain—hips	
Joint pain		Weight gain—waist	
Joints swelling	HHHHH	Wheezing	
Lightheaded when standing			

PLEASE TEAR OFF THIS SHEET AND RETURN IT WITH YOUR TESTING KIT. YOU MUST ALSO RETURN THE BILLING/PHI FORM AND A COPY OF THE FRONT AND BACK OF YOUR INSURANCE CARD TO ENSURE YOUR SAMPLES ARE PROCESSED.

FEMALE: FIRST MORNING VOID COLLECTION INSTRUCTIONS

IMPORTANT: Confirm that it is **YOUR NAME** on all collection tubes. Please call 877-316-8686 if your name is **NOT** listed on the collection tubes.

TAKE ALL MEDICATIONS AS PRESCRIBED! See exceptions*

*If you are taking oral DHEA, oral testosterone, oral progesterone and/or oral estrogen

Take 16 hours before collecting. After collection, resume as prescribed.

TIMING GUIDE - CHOOSING YOUR COLLECTION DATE

① If you experience a menstrual cycle:

• The optimal time to collect is 7 days before your next period.

② If you use patches

- Start collecting 1-2 days after applying patch.
- ① If you administer or receive injections/pellets
- Start collecting at midpoint between injection/pellet inserts.

If none of these apply to you Vou may collect on any day unless your

You may collect on any day <u>unless</u> you are taking hormones. If you are taking hormones, collect on the days you take hormones.

If you use cream/topical hormones, melatonin, hydrocortisone, pregnenolone and/or placental derivatives

Start collection at least 16 hours after the last dose.

① If you use cream/topical hormones on genitals:

- Apply your hormones normally until the day prior to testing. The day prior to testing - apply genital hormone cream to your upper inner arm and rub in well.
- The day of testing wake up and use a "Clean Catch" technique to collect your sample, then apply your cream as you normally do after your morning collection.

DO NOT STOP

TAKING YOUR

HORMONES!

 Clean Catch Technique - Wipe from front to back with a damp cloth. Hold the labia apart while collecting urine. This will help prevent the urine from washing the vaginal walls and contaminating the urine with residual cream.

If you are unsure about the exceptions and/or timing guide above, please contact our staff at 877-316-8686 for assistance before collecting.

Instructions for Collection

The night before testing:

- Limit your caffeine & diuretic intake for the next 24 hours. Stop drinking fluids between 8pm and your collection.
- Completely empty your bladder before bedtime. **Do not collect this urine.**
- Place the ice pack into the freezer.

On the day of collection:

- Upon waking, collect urine into a disposable cup (not provided). Transfer urine from cup to yellow-top tube labeled "waking collection" as indicated here using plastic pipette (included).
- Log your urination time on the page below where indicated. You may also enter your collection time online by logging in to PhysiciansLab.com, or by using our email and text messaging instructions.
- Sample can be kept at room temperature for a maximum of 12 hours, after which it MUST be refrigerated.
- Keep urine samples in your refrigerator until they are ready to be shipped.

Shipping:

Place the urine samples and frozen ice pack into the box immediately before shipping and return using the pre-paid shipping envelope provided.

Tear here and return form below with kit

FEMALE Q	UESTIONNAIRE:	GENERAL INFO	RMATION
· · · · · · · · · · · · · · · · · · ·			

You may skip this page if: You complete this questionnaire online at PhysiciansLab.com or through our text messaging and email instructions.
Check here if completed online:
Name:
Height:/ feet/inches Weight: pounds Currently Pregnant? No Yes Due Date:/ mm/dd/yy
Hysterectomy? No Yes Hysterectomy Year: 4 digit year Ovaries Removed? N/A One Both Ovaries Removed Year: 4 digit yr
Menstrual Cycles: Regular Irregular None First Day of Last Menses://
Collection Day:/ mm/dd/yy

Please log your collection time below.

See instruction page for details on when and how to collect, as well as contact information for our customer support team if you have any questions.



Waking Collection:

Collect your sample in the **yellow-top** tube labeled "waking collection" upon waking.

Waking Collection Time ____:___

Medication Questions

Currently taking inflammatory medication containing cortisone or other anti-inflammatory steroids?

No Yes I am not currently taking any hormones

Hormone Medications

The medication information collected in the section below is essential for the accuracy of your laboratory results and will be displayed therein. **Please check off all applicable information.** If you select one or more Delivery Method(s) please indicate how long you have been taking each hormone.

Delivery Method	Testosterone*	Progesterone*	Estradiol*	DHEA*	Pregnenolone*	Estriol*	Adrenal Supplement*
Oral							
Injection							
Pellet							
Cream/Gel							
Troche							
Sublingual							
Patch							
How Long?							
0-3 Months							
3-6 Months							
6-12 Months							
12 .							

^{*}If you are taking a medication that combines any of the hormones listed above, please check each of those individual hormones.



Continued on back